

Disposable Rf Plasma Electrodes


Indications

This product is used in conjunction with the surgical system. The radio frequency plasma surgical treatment system transmits radio frequency energy to the lesion tissue through this electrode to perform ablation, cutting, coagulation or hemostasis.



Features

- ※ Complete ablation and hemostasis in all directions under transforaminal endoscope, without damaging tissue and nerves
- ※ Controllable power for precise surgery
- ※ Precise ablation of nucleus pulposus tissue within the intervertebral disc
- ※ Short operation time and less pain for the patient
- ※ No sequelae, low recurrence rate



Specifications

Code	Working Length (mm)	Plug Type
RPE08-40S	400 ± 10mm	Type S - Rectangle
RPE08-34S	340 ± 10mm	
RPE08-24S	240 ± 10mm	
RPE08-15S	150 ± 10mm	
RPE08-40R	400 ± 10mm	Type R-Round
RPE08-34R	340 ± 10mm	
RPE08-24R	240 ± 10mm	
RPE08-15R	150 ± 10mm	
RPE08-40F	400 ± 10mm	Type F-Split
RPE08-34F	340 ± 10mm	
RPE08-24F	240 ± 10mm	
RPE08-15F	150 ± 10mm	