

Endoscopic Injection Needle

Indications

This product is used for liquid drug injection in the human digestive tract under endoscope.



Features

- ※ Unique lock reminder
- ※ Ultra-smooth and ultra-permeable catheter
- ※ Easy to observe blood return
- ※ Non-adhesive
- ※ Lightweight and flexible handle design
- ※ Unique stroke control, accurate needle length
- ※ The outer tube sheath is designed with a one-time arc design, which effectively controls the needle tip guidance, and prevents accidents caused by deep needle penetration or follow-up of the sheath during acupuncture operations.



Specifications

Code	Needle Length	Needle Outer Diameter	Sheath Length	Outer Diameter
SN18-07/195	5 mm	0.7mm (22g)	195cm	2.4mm
SN18-05/160	5mm	0.5mm (25g)	160cm	2.4mm
SN18-07/160	5mm	0.7mm (22g)	160cm	2.4mm
SN18-09/160	5mm	0.9mm (20g)	160cm	2.4mm
SN18-05/200	5mm	0.5mm (25g)	200cm	2.4mm
SN18-07/200	5mm	0.7mm (22g)	200cm	2.4mm
SN18-09/200	5mm	0.9mm (20g)	200cm	2.4mm
SN18-05/230	5mm	0.5mm (25g)	230cm	2.4mm
SN18-07/230	5mm	0.7mm (22g)	230cm	2.4mm
SN18-09/230	5mm	0.9mm (20g)	230cm	2.4mm
SN19-05/230	5mm	0.5mm (25g)	230cm	1.8mm
SN19-07/230	5mm	0.7mm (22g)	230cm	1.8mm